

Sample Collection Set



Stable
prevent nucleic acid degradation

Safe
inactivated, reduce infectivity and pathogenicity

Convenient
storage under room temperature and convenient for transportation

User-friendly
with different medium specification to meet diverse needs



The sample collection set is used for storage and transportation of human specimen. It offers biological buffer, maintains normal PH value of organism and stability of nucleic acid and prevents nucleic acid degradation. Besides, the product can inactivate the advanced structure of pathogen protein, deprive biological activity and reduce infectivity and pathogenicity.

ORDERING INFORMATION


Product Name	Sample Collection Set		
Specification	100T/kit		
Product Code	T380H	T381H	T382H
Preservation Solution	3mL	6mL	12mL
	Individual testing	Pool testing (10 in 1)	Pool testing (20 in 1)
Specimen Collection Method	Nasopharyngeal swab and Oropharyngeal swab		
Storage and Stability	Store under room temperature; Validity period:24 months		
Usage	Storage and transportation of human specimen		

TEST PROCEDURE


Oral sampling




① Swab on the posterior wall of the pharynx tonsil recess, lateral wall, etc. Repeatedly wipe 3-5 times.



② Immerse the swab head after sampling preserved in liquid in the tube.



③ Break off the swab from the breaking point (30mm)




④ Tighten the lid and shake

Nasal sampling



① Swab along the nasal passage to the back wall of the nasopharyngeal cavity, gently rotate and slowly remove the swab.



② Immerse the swab head after sampling preserved in liquid in the tube.



③ Break off the swab from the breaking point (80mm/52mm)



④ Tighten the lid and shake

